



Analysis of organochlorine pesticides by Gas Chromatograph Chromatec-Crystal 5000 with MSD

Summary

Organochlorine (OC) pesticides were widely used in agriculture and pesticides control until research and public concern regarding the hazards of their use led to government restrictions and bans. Despite restrictions and bans on the use of many organochlorine pesticides in the 1970s and 1980s, they continue to persist in the environment today.

GC-MSD is the most effective instrument for OCP environment control.

Analysis methods

1. US EPA Method 525. Determination of semivolatile organic chemicals in drinking water by solid phase extraction and capillary column gas chromatography / mass spectrometry (GC/MS).

Instrument configuration

- Gas chromatograph Chromatec-Crystal 9000
- Split-splitless (SSI) inlet or Programmable Split/Splitless Inlet (PSSI)
- Column BP-5ms (30 m × 0.25 mm × 0.25 μm), Cat. # 054310
- Detector – MSD
- Carrier-gas – helium

Analysis mode

Chromatograph

Analysis time 25 min

Column

Carrier-gas flow 1 ml/min

Column temperature

Isotherm 1: 80 °C 1 min 20 °C/min

Isotherm 2: 180 °C 0 min 5 °C/min

Isotherm 3: 270 °C

Split/splitless Inlet

Temperature 250 °C

Injection Mode Splitless

MSD

Ion source temperature 300 °C

Transfer line temperature 300 °C

Scan Mode Full Scan

Mass Range 45 - 550 amu


Results

Figure 1 shows the chromatogram of organochlorine pesticides in TIC mode.

The system was calibrated at 6 levels: 20, 50, 100, 200, 500, 100 ppb. Each calibration level contained 26 pesticides and internal standard (anthracene-d10). The calibration was performed in SIM mode. The chromatogram in SIM mode for 50 ppb level is shown in Figure 2.

Figure 3 shows the calibration curves of some components.

Chromatograms



- | | | |
|----------------------------|--------------------------|------------------------|
| 1. α -BHC | 10. Heptachlor epoxide A | 19. Endrin |
| 2. Hexachlorobenzene | 11. Heptachlor epoxide B | 20. Endosulfan II |
| 3. β -BHC | 12. cis-Chlordane | 21. o,p'-DDT |
| 4. γ -BHC | 13. o,p'-DDE 1 | 22. p,p'-DDD |
| 5. Pentachloronitrobenzene | 14. Endosulfan I | 23. Endrin aldehyde |
| 6. Anthracene-d10 | 15. trans-Chlordane | 24. Endosulfan sulfate |
| 7. δ -BHC | 16. p,p'-DDE | 25. p,p'-DDT |
| 8. Heptachlor | 17. Dieldrin | 26. Endrin ketone |
| 9. Aldrin | 18. o,p'-DDD | 27. Methoxychlor |

Figure 1 – Organochlorine pesticides mix (500 ppb) using Scan Mode




Figure 2 – Organochlorine pesticides mix (50 ppb) using SIM Mode




Figure 3 – Calibration curves of some pesticides

Conclusion

Analysis of organochlorine pesticides can be accomplished using GC-MS "Chromatec-Crystal 9000". To achieve high sensitivity, the analysis is performed in SIM mode. The calibration range is from 20 to 1000 ppb.